


WILAméd

Equipment for Professionals



NEO Hi-Flow Nasal Cannula

www.wilamed.com



NEO Hi-Flow Nasal Cannula (NHFT)

The use of **Hi-Flow Nasal Cannula** is increasingly popular for non-invasive respiratory support in Neonatal Intensive Care Units (NICU's) around the world.

The specially designed WILamed NEO Hi-Flow Nasal Cannula for nasal Hi-Flow O₂ Therapy (NHFT) can deliver heated and humidified gas at flow rates of up to 8 liter per minute .

They are a simpler, lighter and a more comfortable interface with the infant and advantageous for Mother-Infant bonding (Kangaroo Care).

Order Details

Order number	Item description	Max. Flow Rate	QTV.
270.995	PREMATURE Hi-Flow Nasal Cannula	6 LPM	25/Box
270.996	NEONATE Hi-Flow Nasal Cannula	7 LPM	25/Box
270.997	INFANT Hi-Flow Nasal Cannula	7 LPM	25/Box
270.998	PEDIATRIC Hi-Flow Nasal Cannula	8 LPM	25/Box



Distributed by:

WILamed GmbH

Medizinische Geräte und Zubehör

Gewerbepark Barthelmesaurach

Aurachhöhe 5-7

91126 Kammerstein (Germany)

Phone: +49 9178 996999-0

Fax: +49 9178 996778

info@wilamed.com

www.wilamed.com

NEO Hi-Flow